IMPORTANT INFORMATION – PLEASE READ BEFORE USE  
INFORMATIONS IMPORTANTES – À LIRE AVANT UTILISATION  
INFORMAZIONI IMPORTANTI – LEGGERE PRIMA DELL'USO  
WICHTIGE INFORMATIONEN – VOR GEBRAUCH BITTE LESEN  
INFORMACIÓN IMPORTANTE – LÉALA ANTES DE UTILIZAR ESTE PRODUCTO  
重要な情報 – 使用前にお読みください  
ВАЖНА ИНФОРМАЦИЯ – МОЛЯ, ПРОЧЕТЕТЕ ПРЕДИ УПОТРЕБА  
VAŽNE INFORMACIJE – PROČITATI PRIJE UPORABE  
DŮLEŽITÉ INFORMACE – PŘEČTĚTE SI PŘED POUŽITÍM  
VIGTIGE OPLYSNINGER – SKAL LÆSES INDEN BRUG  
BELANGRIJKE INFORMATIE – VOOR GEBRUIK DOORLEZEN  
OLULINE TEAVE – PALUN LUGEGE ENNE KASUTAMIST  
TÄRKEITÄ TIETOJA – LUETTAVA ENNEN KÄYTTÖÄ  
ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ – ΠΑΡΑΚΑΛΟΥΜΕ ΔΙΑΒΑΣΤΕ ΠΡΙΝ ΑΠΟ ΤΗ ΧΡΗΣΗ  
FONTOS TUDNIVALÓK – HASZNÁLAT ELŐTT OLVASSA EL  
SVARĪGA INFORMĀCIJA – LŪDZU, IZLASIET PIRMS LIETOŠANAS  
SVARBI INFORMACIJA – PERSKAITYKITE PRIEŠ NAUDODAMI  
VIKTIG INFORMASJON – LES FØR BRUK  
WAŽNE INFORMACJE – NALEŻY PRZECZYTAĆ PRZED UŻYCIEM  
INFORMAÇÕES IMPORTANTES – LER ANTES DE UTILIZAR  
INFORMAȚII IMPORTANTE – VĂ RUGĂM SĂ LE CITIȚI ÎNAINTEA UTILIZĂRII  
DŮLEŽITÉ INFORMÁCIE – PREČÍTAJTE SI ICH PŘED POUŽITÍM  
ПОМЕМБНЕ ИНФОРМАЦИЈЕ – ПРЕБЕРИТЕ ПРЕД УПОРАБО  
VIKTIG INFORMATION – LÄS FÖRE ANVÄNDNING  
ÖNEMLI BILGI – LÜTFEN KULLANMADAN ÖNCE OKUYUN

Contents (English - EN).....	1
Contenu (Français - FR).....	7
Contenuto (Italiano - IT) .....	14
Inhalt (Deutsch - DE) .....	21
Índice (Español - ES).....	28
内容 (日本語 - JA).....	35
Съдържание (Български - BG).....	41
Sadržaj (Hrvatski - HR).....	48
Obsah (Česky - CS).....	54
Indhold (Dansk - DA).....	60
Inhoud (Nederlands - NL).....	66
Sisukord (Eesti - ET).....	73
Sisältö (Suomi - FI) .....	79
Περιεχόμενα (Ελληνικά - EL).....	85
Tartalom (Magyar - HU) .....	92
Saturs (Latviski - LV).....	99
Turinys (Lietuvių - LT).....	105
Innhold (Norsk - NO).....	111
Spis treści (Polski - PL) .....	117
Índice (Português - PT).....	124
Cuprins (Română - RO).....	131
Obsah (Slovenčina - SK).....	138
Vsebina (Slovenščina - SL).....	144
Innehåll (Svenska - SV) .....	150
İçindekiler (Türkçe - TR).....	157

## **Instructions for Use/Intended Use**

Tecomet delivery systems (cases and trays) are intended for the protection, organization, and delivery to the surgical field of surgical instruments and/or other medical devices. Tecomet delivery systems (cases and trays) are not designed to maintain sterility by themselves. They are designed to facilitate the sterilization process when used in conjunction with a wrapping material (FDA cleared sterilization wrap). Wrapping material is designed to allow air removal, steam penetration/evacuation (drying) and maintain the sterility of the internal components.

## **Warnings and Precautions**

### **Warnings**

- The reusable instruments are provided NON-STERILE  and must be properly cleaned and sterilized prior to each use.
- The delivery system is not designed to maintain sterility by itself. Only use the delivery system with an approved sterilization wrapping material. The container has feet that are designed for storage purposes only and not to stack cases for sterilization.
- Read these instructions completely before using the devices.
- It is the surgeon's responsibility to be familiar with the appropriate surgical techniques prior to use of the devices.
- Personal Protective Equipment (PPE) should be worn when handling or working with contaminated or potentially contaminated instruments.
- Do not use metal brushes, abrasive cleansers or abrasive pads for cleaning.

### **Precautions**

- U.S Federal law restricts this device to sale by or on the order of a physician.
- The instruments should be inspected for damage and dull/uneven cutting edges prior to each use. Instruments that show signs of damage or degradation should not be used.

- The sterilization container should be inspected for damage prior to each use. Careful attention should be paid to the handle integrity and the container exterior to ensure no sharp edges/burrs are present to tear the wrap.

### **Single Use Instruments**

- Risks of reusing single use instruments include but are not limited to infection to the patient and/or reduced reliability of function.

### **Limitations of Reprocessing**

Repeated processing according to these instructions has minimal effect on the instruments. End of life for stainless steel or other metal surgical instruments is normally determined by wear and damage due to the intended use or misuse and not reprocessing.

### **Materials and Restricted Substance**

- For indication that the device contains a restricted substance or material of animal origin see product label.

### **Cleaning**

- It is highly recommended that as soon as possible after each surgical procedure, and before sterilization, the instruments be cleaned in order to limit the drying time of residue biologic soil left on the instruments.
- Water quality used for diluting cleaning agents and for rinsing instruments should be carefully considered. Use of distilled water for cleaning and sterile water for rinsing is recommended. Avoid using hot water as this will coagulate and harden protein based soil.
- All cleaning agents and disinfectants must be prepared according to the recommendations of their manufacturer. Only use cleaning agents and disinfectants that have a nearly neutral pH and are approved for use on surgical instruments.

### **Point of Use Pre-cleaning**

- Remove excess biologic soil and tissue from instruments using disposable wipes.
- As soon as possible after use, set instruments in a basin of distilled water or in a tray covered with damp towels.

## **Manual Cleaning of Instruments**

1. Prepare a proteolytic enzyme based detergent like Enzol (or equivalent) and prepare a cleaning solution according to the detergent manufacturer's recommendation.
2. Immerse instruments and soak for the time recommended by the detergent manufacturer.
  - While immersed actuate all hinged devices and those with moving parts.
3. Use a soft bristle cleaning brush and scrub the instruments until all visible contamination has been removed. Scrub the device below the surface of the cleaning solution to prevent aerosolization of contaminants. Pay particular attention to the features that will pose a challenge to effective cleaning. Follow these general guidelines:
  - Using the brush, scrub all surfaces of the device. Pay particular attention to any cutting features or rough surfaces used for filing or abrading.
  - Use a snug fitting pipe cleaner to access blind holes or cannula. Insert and rotate in blind holes or pass the pipe cleaner through any cannula three (3) times.
4. Rinse thoroughly with distilled water until all traces of cleaning solution are removed.
  - Actuate all hinged devices and those with moving parts while rinsing.
5. Prepare an ultrasonic bath with a cleaning solution at the concentration and temperature recommended by the detergent manufacturer.
6. Submerge the instruments and activate the bath for a minimum of 10 minutes. A frequency of 25 – 50 kHz is recommended.
7. Remove and rinse the instruments in distilled or sterile water for at least one (1) minute or until all traces of cleaning solution are removed.
  - Actuate all hinged devices and those with moving parts while rinsing.
8. Visually inspect the instruments for visible soil and repeat these cleaning steps if remaining soil is observed.
9. Dry the instruments with clean, lint-free wipes in preparation for sterilization. Use clean pressurized air to remove moisture from hard to reach areas.

## **Automated Cleaning of Instruments Using Washer-Disinfector**

1. Prepare a solution of enzymatic detergent according to the manufacturer's recommendations.
2. Immerse instruments and soak for the time recommended by the detergent manufacturer.
  - While immersed actuate all hinged devices and those with moving parts.
3. Use a soft bristle cleaning brush and scrub the instruments until all visible contamination has been removed. Scrub the device below the surface of the cleaning solution to prevent aerosolization of contaminants. Pay particular attention to the features that will pose a challenge to effective cleaning. Follow these general guidelines:
  - Using the brush scrub all surfaces of the device. Pay particular attention to any cutting features or rough surfaces used for filing or abrading.
  - Use a snug fitting pipe cleaner to access blind holes or cannula. Insert and rotate in blind holes or pass the pipe cleaner through any cannula three (3) times.
4. Rinse thoroughly with distilled water until all traces of cleaning solution are removed.
  - Actuate all hinged devices and those with moving parts while rinsing.
5. Load instruments in an automated washer-disinfector in a manner that maximizes exposure of the instrument surfaces.
6. Operate the washer-disinfector according to the manufacturer's instructions to ensure all cycle parameters (i.e. time, temperature) are followed.
7. Remove instruments and check for remaining soil or wetness. If soil remnants are observed repeat the automated cleaning cycle. If remaining wetness is observed dry the instruments with clean, lint-free wipes in preparation for sterilization.

## **Manual Cleaning of Sterilization Container**

1. Prepare a solution of neutral pH detergent according to the manufacturer's recommendations. Cleansers that are too acidic or alkaline may damage the anodized finish of the aluminum case.
2. Using a soft sponge or cloth, clean all surfaces of the container base, lid and instrument trays.

3. Thoroughly rinse the container components under clean running water to remove all residue detergent.
4. Thoroughly dry the container components.

### **Automated Cleaning of Sterilization Container Using Washer-Disinfector**

1. Prepare a solution of neutral pH detergent according to the washer manufacturer's recommendations.
2. Place the container components into the washer in a manner that will prevent them from moving and start the cycle.
3. After the cleaning cycle is complete remove the container components and verify they are dry. If wetness is observed, dry the components with clean, lint-free wipes.

### **Sterilization**

Moist heat/steam sterilization is the preferred and recommended method for the Tecomet delivery systems.

The instruments and container must be properly cleaned before sterilization.

Place the instruments in their respective position within the sterilization container according to the markings / labeling in the container. Once the container is loaded put the lid on and secure all lid locks. The closed container should now be wrapped with an approved wrap according to the AAMI ST79 recommendations. The container is now ready for sterilization according to the parameters listed below.

#### Recommended Parameters for Steam Sterilization

<b>Sterilization Mode</b>	<b>Exposure Temp.</b>	<b>Exposure Time</b>	<b>*Dry Time</b>
Dynamic-Air-Removal	132°C (270°F)	4 minutes	60 minutes
Dynamic-Air-Removal	135°C (275°F)	3 minutes	60 minutes

\* Fully loaded cases validated with 60 minutes cooling time using two applications of Kimguard KC600 single ply wrap using sequential envelope technique per AAMI ST79. The drying time may vary due to sterilization equipment used, wrapping method, and material. It is the health care facility's responsibility to validate the appropriate drying time using their process.

## Labeling Symbols

Symbol for “Manufacturer”

Catalog Number

Caution. Consult Accompanying Documents

Consult Instructions for Use

Lot Number

Non-Sterile Instrument

U.S Federal law restricts this device to sale by or on the order of a physician

Sterile Instrument

CE <sup>1</sup>

CE Mark with Notified Body Number<sup>1</sup>

Authorized Representative in the European Community

Medical Device

Packaging Unit

Country of Manufacture

Distributor

<sup>1</sup>Refer to the labeling for CE information

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Symmetry Medical Manufacturing, Inc.
Symmetry Medical Polyvac S.A.S