IMPORTANT INFORMATION - PLEASE READ BEFORE USE
INFORMATIONS IMPORTANTES - À LIRE AVANT UTILISATION
INFORMAZIONI IMPORTANTI - LEGGERE PRIMA DELL'USO
WICHTIGE INFORMATIONEN - VOR GEBRAUCH BITTE LESEN
INFORMACIÓN IMPORTANTE - LÉALA ANTES DE UTILIZAR ESTE PRODUCTO
重要な情報 - 使用前にお読みください
ВАЖНА ИНФОРМАЦИЯ - МОЛЯ, ПРОЧЕТЕТЕ ПРЕДИ УПОТРЕБА
VAŽNE INFORMACIJE - PROČITATI PRIJE UPORABE
DŮLEŽITÉ INFORMACE - PŘEČTĚTE SI PŘED POUŽITÍM
VIGTIGE OPLYSNINGER - SKAL LÆSES INDEN BRUG
BELANGRIJKE INFORMATIE - VOOR GEBRUIK DOORLEZEN
OLULINE TEAVE - PALUN LUGEGE ENNE KASUTAMIST
TÄRKEITÄ TIETOJA - LUETTAVA ENNEN KÄYTTÖÄ
ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ - ΠΑΡΑΚΑΛΟΥΜΕ ΔΙΑΒΑΣΤΕ ΠΡΙΝ ΑΠΟ
ΤΗ ΧΡΗΣΗ
FONTOS TUDNIVALÓK - HASZNÁLAT ELŐTT OLVASSA EL
SVARĪGA INFORMĀCIJA - LŪDZU, IZLASIET PIRMS LIETOŠANAS
SVARBI INFORMACIJA - PERSKAITYKITE PRIEŠ NAUDODAMI
VIKTIG INFORMASJON - LES FØR BRUK
WAŽNE INFORMACJE - NALEŻY PRZECZYTAĆ PRZED UŻYCIEM
INFORMAÇÕES IMPORTANTES - LER ANTES DE UTILIZAR
INFORMAȚII IMPORTANTE - VĂ RUGĂM SĂ LE CITIȚI ÎNAINTEA UTILIZĂRII
ВАЖНАЯ ИНФОРМАЦИЯ - ИЗУЧИТЕ ПЕРЕД ПРИМЕНЕНИЕМ ИЗДЕЛИЯ
VAŽNE INFORMACIJE - PROČITATI PRIJE UPORABE
DŮLEŽITÉ INFORMÁCIE - PREČÍTAJTE SI ICH PRED POUŽITÍM
ROMEMBNE INFORMACIJE - PREBERITE PRED UPORABO
VIKTIG INFORMATION - LÄS FÖRE ANVÄNDNING
ÖNEMLİ BİLGİ - LÜTFEN KULLANMADAN ÖNCE OKUYUN

Contents (English - EN).....	1
Contenu (Français - FR).....	5
Contenuto (Italiano - IT).....	9
Inhalt (Deutsch - DE).....	13
Índice (Español - ES).....	18
内容 (日本語 - JA).....	22
Съдържание (Български - BG).....	26
Sadržaj (Hrvatski - HR).....	31
Obsah (Česky - CS).....	35
Indhold (Dansk - DA).....	39
Inhoud (Nederlands - NL).....	43
Sisukord (Eesti - ET).....	47
Sisältö (Suomi - FI).....	51
Περιεχόμενα (Ελληνικά - EL).....	55
Tartalom (Magyar - HU).....	60
Saturs (Latviski - LV).....	64
Turinys (Lietuvių - LT).....	68
Innhold (Norsk - NO).....	72
Spis treści (Polski - PL).....	76
Índice (Português - PT).....	80
Cuprins (Română - RO).....	84
Содержание (Русский - RU).....	88
Sadržaj (Srpski - SR).....	93
Obsah (Slovenčina - SK).....	97
Vsebina (Slovenščina - SL).....	101
Innehåll (Svenska - SV).....	105
İçindekiler (Türkçe - TR).....	109

Instructions for Use/Intended Use

The X Series Power System Instrument Cases are intended for the protection, organization, and delivery to the surgical field of surgical instruments and/or other medical devices. The X Series Power System Cases are not intended to maintain sterility by itself. It is designed to facilitate the sterilization process when used with an FDA approved sterilization wrap or rigid sterilization container. Wrapping material is designed to allow air removal, steam penetration/evacuation (drying) and to maintain the sterility of the internal components.

Indications For Use

The X-Series Cases are indicated for use to protect, organize, transport, and facilitate the sterilization process by allowing steam penetration and air removal. When used in conjunction with an approved sterilization wrap and following the specified IFU, the sterility of the enclosed medical device is maintained provided the sterilized instruments are properly housed and sterility maintenance is performed per *ANSI ST79*.

Intended Patient Population

None. The cases are intended to protect, organize, transport, and facilitate the sterilization the X-Series Power System Instruments.

Intended User

The X-Series Cases are intended to be used by qualified health care professionals e.g., surgeons, physicians, operating room staff, and individuals involved in preparation of the surgical device.

Expected Clinical Benefits

The cases protect, organize, transport, and facilitate the sterilization the X-Series Power System Instruments.

Performance Characteristics

The cases protect, organize, transport, and facilitate the sterilization the X-Series Power System Instruments.

Adverse Events & Complications

All surgical operations carry risk. The following are frequently encountered adverse events and complications related to having a surgical procedure in general:

- Delay to surgery caused by missing, damaged or worn instruments.
- Tissue injury and additional bone removal due to blunt, damaged or incorrectly positioned instruments.
- Infection and toxicity due to improper processing.

Adverse events to user:

- Cuts, abrasions, contusions or other tissue injury caused by burs, sharp edges, impaction, vibration or jamming of instruments.

Adverse Events & Complications – Reporting of Serious Incidents

Serious Incident Reporting (EU)

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- The death of a patient, user, or other person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

Materials & Restricted Substances

For indication that the device contains a restricted substance or material of animal origin see product label.

Combination of Medical Devices

The cases are designed to protect, organize, transport, and facilitate the sterilization the X-Series Power System Instruments. The cases are not to be used in combination with devices provided by other manufacturers. Tecomet assumes no liability in such instances. Make sure all latches and handles are secure and in working order, and the instruments are placed in their intended locations/brackets.

Recommendations for Care, Cleaning, and Sterilization for the X Series Power System Instruments Case

Both physical and chemical (detergent) processes may be necessary to clean soiled items. Chemical (detergent) cleaners alone cannot remove all soil and debris; therefore, a careful manual cleaning of each item with a soft sponge or cloth is essential for maximum decontamination. For difficult areas, a clean soft bristled brush is recommended. Once the items have been cleaned, they should be thoroughly rinsed with clean water to remove any detergent or chemical residue before sterilization. **Tecomet** recommends the use of a mild enzymatic detergent with a near neutral pH. Do not use solvents, abrasive cleaners metal brushes, or abrasive pads. Baskets may be placed in mechanical cleaning equipment. **Always** inspect instruments and case for cleanliness and damage before use. Make sure all latches and handles are secure and in working order. Do not overload baskets. Tecomet recommends that the baskets be wrapped according to the sterilization wrap manufacturer's instructions prior to sterilization to maintain sterility of internal components/items and for proper aseptic presentation to the surgical field. Always follow the X Series Owner's Manual if their sterilization or drying recommendations exceed these guidelines. For steam autoclaves, Tecomet has verified product performance in the following cycles.

Pre-vacuum Sterilizer

Wrapped baskets should be sterilized according to the table below.

Cycle Type	Temperature	Exposure Time	Dry Time	Cool-Down Time*
United States Recommended Parameters				
Pre-vacuum / Vacuum Pulse	132°C / 270°F	4 minutes	45 minutes	30 minutes
Cycle Type	Temperature	Exposure Time	Dry Time	Cool-Down Time*
European Recommended Parameters				
Pre-vacuum/ Vacuum Pulse	134°C / 273°F	3 minutes	45 minutes	30 minutes

*Allow the wrapped device to cool outside of the autoclave for 30 minutes prior to use.

Filtered Sterilization Container Systems

Tecomet delivery systems have not been validated for use in filtered sterilization containers. Tecomet does not recommend use of such a system. If filtered sterilization container systems are used, the user is responsible for following the recommendations of the manufacturer for the proper placement and use of cases and tray inside the container.

Contraindications

The X Series Power System Instrument Cases have **NOT** been validated for use with flexible endoscopes, or devices with lumens or working channels longer than 4 inches (3mm ID). Always refer to the X Series Owner's Manual. The X Series Power System

Instrument Cases have **NOT** been validated for ETO sterilization of devices. Always refer to the X Series Owner's Manual. Stacking the X Series Power System Instrument Cases and overloading the cases will adversely affect sterilization and drying effectiveness. **DO NOT STACK** cases or trays in an autoclave chamber.

Notice

1. **DO NOT** load cases into sterilizer on sides or upside down with lid side on shelf or cart. Load cases on cart or shelf so that the lid is always facing upwards. This will allow for proper drying. The X Series Power System Instrument Cases are designed to drain in this position.
2. **DO NOT STACK** cases or trays in an autoclave chamber.
3. After the autoclave door is opened, all cases must be allowed to cool thoroughly. Place cases on rack or shelf with linen cover until cooling is complete. The potential for condensation may increase if the case is not allowed to cool properly.
4. If condensation is observed, check to ensure steps 1, 2, and 3 were followed. In addition, verify that the steam, which is being used for sterilization processing, has a quality of more than 97%¹. Also confirm that the sterilizers have been inspected for routine maintenance in accordance with manufacturer recommendations.

Disposal

At the end of the device's life safely dispose of the device in accordance with local procedures and guidelines. Any device that has been contaminated with potentially infectious substances of human origin (such as bodily fluids) should be handled according to hospital protocol for infectious medical waste. Any device that contains sharp edges should be discarded according to hospital protocol in the appropriate sharps container.

SYMBOLS USED ON LABELING

Caution

Non-Sterile

U.S Federal law restricts this device to sale by or on the order of a physician

CE Mark

Authorized Representative in the European Community

Manufacturer

Date of Manufacture

Lot Number

Catalog Number

Consult Instructions for Use

European Union Importer

Medical Device

Unique Device Identifier

Packaging Unit

Country of Manufacture

Distributor

Swiss Authorized Representative

The following table provides guidance on maximum loads. Do not overload delivery systems (baskets), and always follow AAMI, OSHA and hospital standards for maximum loading. Follow the X Series Owner's Manual's instructions for use.

Tray Size	Maximum Total Weight
Small (~ 4 x 8 x 1 in.)	1 lb. (0.45 kg)
Tray (~ 8 x 17 x 1 in.)	2 lb. (0.91 kg)
½ Size Case (~ 7 x 7 x 4 in.)	8 lb. (3.6 kg)
¾ Size Case (~ 7 x 12 x 4 in.)	13 lb. (5.9 kg)
Full Size Case (~ 9 x 19 x 4 in.)	22 lb. (10.0 kg)

REFERENCES

¹ANSI/AAMI ST79: Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities

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Symmetry Medical Manufacturing, Inc.
Symmetry Medical Polyvac S.A.S