IMPORTANT INFORMATION – PLEASE READ BEFORE USE
INFORMATIONS IMPORTANTES – À LIRE AVANT UTILISATION
INFORMAZIONI IMPORTANTI – LEGGERE PRIMA DELL'USO
WICHTIGE INFORMATIONEN – VOR GEBRAUCH BITTE LESEN
INFORMACIÓN IMPORTANTE – LÉALA ANTES DE UTILIZAR ESTE PRODUCTO

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Shoulder Revision Extraction Instruments

INSTRUCTIONS FOR USE

Description

The Shoulder Revision Extraction Instrument set is comprised of a complete set of hand held instruments required for the revision of a previous shoulder arthroplasty. The instruments are selected and used at the discretion of the surgeon to perform various surgical actions dependent upon the specific device function. The set is provided with single-use pre-sterile osteotomes and a sterilization container that holds the remaining reusable instruments comprised of chisels, pliers, cement separators and splitters, curettes, rongeurs, tamps, mallets, clamps, extractors, slap-hammers, drill guides and bits, lamina spreaders, forks, and handles.

Intended Use

The Shoulder Revision Extraction Instruments made by Tecomet includes single-use sterile osteotome blades and a set of reusable hand held instruments (e.g. curettes, pliers, chisels, extractors, cement separators, clamps, mallet, etc.) designed to provide various capabilities for removing implants from a previous shoulder arthroplasty. The system also includes a sterilization container designed to organize, protect and facilitate sterilization of the reusable instruments when used in conjunction with an approved wrapping material (e.g. FDA cleared sterilization wrap).

Intended Patient Population

The device is prescriptive, therefore, a knowledgeable orthopedic surgeon may utilize the device on any patient they deem necessary. The device is to be used on patients undergoing total shoulder arthroplasty revision.

Indications For Use

The Shoulder Revision Extraction Instrument set is comprised of a complete set of hand held instruments required for the revision of a previous shoulder arthroplasty. The instruments are selected and used at the discretion of the surgeon to perform various surgical actions dependent upon the specific device function. The set is provided with single-use pre sterile osteotomes and a sterilization container that holds the remaining reusable instruments comprised of chisels, pliers, cement separators and splitters, curettes, rongeurs, tamps, mallets, clamps, extractors, slap-hammers, drill guides and bits, lamina spreaders, forks, and handles.

Contra-indications

These surgical instruments are prescription use. The instruments are only to be used by qualified health care personnel. There are no other contra-indications for the Shoulder Revision Extraction Instruments.

Intended User

The Shoulder Revision Extraction Instrument set prescriptive and therefore to be used by qualified orthopedic surgeons trained in the respective surgical technique.

Expected Clinical Benefits

When used as intended, Shoulder Revision Extraction Instrument set aids in the safe removal of the previously failed implant.

Performance Characteristics

The Shoulder Revision Extraction Instrument set aids in the safe removal of the previously failed implant.

Adverse Events & Complications

All surgical operations carry risk. The following are frequently encountered adverse events and complications related to having a surgical procedure in general:

- Delay to surgery caused by missing, damaged or worn instruments.
- Tissue injury and additional bone removal due to blunt, damaged or incorrectly positioned instruments.
- Infection and toxicity due to improper processing.

Adverse events to user:

- Cuts, abrasions, contusions or other tissue injury caused by burs, sharp edges, impaction, vibration or jamming of instruments.

Adverse Events & Complications – Reporting of Serious Incidents

Serious Incident Reporting (EU)

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- The death of a patient, user, or other person,
- The temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- A serious public health threat.

WARNINGS AND PRECAUTIONS

Warnings

- The reusable instruments are provided NON-STERILE  and must be properly cleaned and sterilized prior to each use.
- The sterilization container is not designed to maintain sterility by itself. Only use the container system with an approved sterilization wrapping material. The container has feet that are designed for storage purposes only and not to stack cases for sterilization.
- Osteotome blades are provided single-use sterile and should not be re-used.
- Read these instructions completely before using the devices.
- It is the surgeon's responsibility to be familiar with the appropriate surgical techniques prior to use of the devices.
- Personal Protective Equipment (PPE) should be worn when handling or working with contaminated or potentially contaminated instruments.
- Do not use metal brushes, abrasive cleansers or abrasive pads for cleaning.
- Ensure firm connection with the assembled device prior to use.

Precautions

- U.S Federal law restricts this device to sale by or on the order of a physician.
- The instruments should be inspected for damage and dull/uneven cutting edges prior to each use. Instruments that show signs of damage or degradation should not be used.
- The sterilization container should be inspected for damage prior to each use. Careful attention should be paid to the handle integrity and the container exterior to ensure no sharp edges/burrs are present to tear the wrap.

Materials and Restricted Substance

- For indication that the device contains a restricted substance or material of animal origin see product label.

Single Use Instruments

- The osteotome blades are labeled **Single Use** and are intended to be used one time then discarded.
- Risks of reusing single use instruments include but are not limited to infection to the patient and/or reduced reliability of function.

Disposal

- At the end of the device's life safely dispose of the device in accordance with local procedures and guidelines.
- Any device that has been contaminated with potentially infectious substances of human origin (such as bodily fluids) should be handled according to hospital protocol for infectious medical waste. Any device that contains sharp edges should be discarded according to hospital protocol in the appropriate sharps container.

Limitations of Reprocessing

Repeated processing according to these instructions has minimal affect on the instruments. End of life for stainless steel or other metal surgical instruments is normally determined by wear and damage due to the intended use or misuse and not reprocessing.

CLEANING

- It is highly recommended that as soon as possible after each surgical procedure, and before sterilization, the instruments be cleaned in order to limit the drying time of residue biologic soil left on the instruments.
- Water quality used for diluting cleaning agents and for rinsing instruments should be carefully considered. Use of distilled water for cleaning and sterile water for rinsing is recommended. Avoid using hot water as this will coagulate and harden protein based soil.
- All cleaning agents and disinfectants must be prepared according to the recommendations of their manufacturer. Only use cleaning agents and disinfectants that have a nearly neutral pH and are approved for use on surgical instruments.

Point of Use Pre-cleaning

- Remove excess biologic soil and tissue from instruments using disposable wipes.
- As soon as possible after use, set instruments in a basin of distilled water or in a tray covered with damp towels.

Manual Cleaning of Instruments

1. Prepare a proteolytic enzyme based detergent like Enzol (or equivalent) and prepare a cleaning solution according to the detergent manufacturer's recommendation.
2. Immerse instruments and soak for the time recommended by the detergent manufacturer.
 - While immersed actuate all hinged devices and those with moving parts.
3. Use a soft bristle cleaning brush and scrub the instruments until all visible contamination has been removed. Scrub the device below the surface of the cleaning solution to prevent aerosolization of contaminants. Pay particular attention to the features that will pose a challenge to effective cleaning. Follow these general guidelines:
 - Using the brush, scrub all surfaces of the device. Pay particular attention to any cutting features or rough surfaces used for filing or abrading.
 - Use a snug fitting pipe cleaner to access blind holes or cannula. Insert and rotate in blind holes or pass the pipe cleaner through any cannula three (3) times.

4. Rinse thoroughly with distilled water until all traces of cleaning solution are removed.
 - Actuate all hinged devices and those with moving parts while rinsing.
5. Prepare an ultrasonic bath with a cleaning solution at the concentration and temperature recommended by the detergent manufacturer.
6. Submerge the instruments and activate the bath for a minimum of 10 minutes. A frequency of 25 – 50 kHz is recommended.
7. Remove and rinse the instruments in distilled or sterile water for at least one (1) minute or until all traces of cleaning solution are removed.
 - Actuate all hinged devices and those with moving parts while rinsing.
8. Visually inspect the instruments for visible soil and repeat these cleaning steps if remaining soil is observed.
9. Dry the instruments with clean, lint-free wipes in preparation for sterilization. Use clean pressureized air to remove moisture from hard to reach areas.

Automated Cleaning of Instruments Using Washer-Disinfector

1. Prepare a solution of enzymatic detergent according to the manufacturer's recommendations.
2. Immerse instruments and soak for the time recommended by the detergent manufacturer.
 - While immersed actuate all hinged devices and those with moving parts.
3. Use a soft bristle cleaning brush and scrub the instruments until all visible contamination has been removed. Scrub the device below the surface of the cleaning solution to prevent aerosolization of contaminants. Pay particular attention to the features that will pose a challenge to effective cleaning. Follow these general guidelines:
 - Using the brush scrub all surfaces of the device. Pay particular attention to any cutting features or rough surfaces used for filing or abrading.
 - Use a snug fitting pipe cleaner to access blind holes or cannula. Insert and rotate in blind holes or pass the pipe cleaner through any cannula three (3) times.
4. Rinse thoroughly with distilled water until all traces of cleaning solution are removed.
 - Actuate all hinged devices and those with moving parts while rinsing.
5. Load instruments in an automated washer-disinfector in a manner that maximizes exposure of the instrument surfaces.
6. Operate the washer-disinfector according to the manufacturer's instructions to ensure all cycle parameters (i.e. time, temperature) are followed.
7. Remove instruments and check for remaining soil or wetness. If soil remnants are observed repeat the automated cleaning cycle. If remaining wetness is observed dry the instruments with clean, lint-free wipes in preparation for sterilization.

Manual Cleaning of Sterilization Container

1. Prepare a solution of neutral pH detergent according to the manufacturer's recommendations. Cleansers that are too acidic or alkaline may damage the anodized finish of the aluminum case.
2. Using a soft sponge or cloth, clean all surfaces of the container base, lid and instrument trays.
3. Thoroughly rinse the container components under clean running water to remove all residue detergent.
4. Thoroughly dry the container components.

Automated Cleaning of Sterilization Container Using Washer-Disinfector

1. Prepare a solution of neutral pH detergent according to the washer manufacturer's recommendations.
2. Place the container components into the washer in a manner that will prevent them from moving and start the cycle.
3. After the cleaning cycle is complete remove the container components and verify they are dry. If wetness is observed, dry the components with clean, lint-free wipes.

STERILIZATION

Moist heat/steam sterilization is the preferred and recommended method for the shoulder revision extraction instrument set.

The instruments and container must be properly cleaned before sterilization.

Place the instruments in their respective position within the sterilization container according to the markings / labeling in the container. Once the container is loaded put the lid on and secure all lid locks. The closed container should now be wrapped with an approved wrap according to the AAMI ST79 recommendations. The container is now ready for sterilization according to the parameters listed below.

Recommended Parameters for Steam Sterilization

Sterilization Mode	Exposure Temp.	Exposure Time	* Dry Time
Dynamic-Air-Removal	132°C (270°F)	4 minutes	30 minutes
Dynamic-Air-Removal	135°C (275°F)	3 minutes	30 minutes

* Fully loaded cases validated with 60 minutes cooling time using two applications of Kimguard KC600 single ply wrap using sequential envelope technique per AAMI ST79. The drying time may vary due to sterilization equipment used, wrapping method, and material. It is the health care facility's responsibility to validate the appropriate drying time using their process.

LABELING SYMBOLS

Manufacturer

Date of Manufacture

Catalog Number

Caution

Consult Instructions for Use

Lot Number

Non-Sterile Instrument

U.S Federal law restricts this device to sale by or on the order of a physician

Sterile Instrument

CE Mark¹

CE Mark with Notified Body Number¹

Authorized Representative in the European Community

Symbol for "do not re-use," "single use," or "use only once"

Medical Device

Packaging Unit

Country of Manufacture

Distributor


Sterilized using irradiation

Do not use if package is damaged and consult instructions for use

Use by

Swiss Authorized Representative²

Importer

Unique Device Identifier

¹Refer to the labeling for CE information

²Refer to the labeling for Swiss Authorized Representative

Additional Labeling:

"MANUAL USE ONLY" Device should not be attached to a power source and is only intended to be handled manually.

"REMOVE FOR CLEANING & STERILIZATION" Device should be disassembled prior to cleaning and sterilization.

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Symmetry Medical Manufacturing Inc.
DePuy Orthopaedics, Inc
Symmetry Medical Polyvac S.A.S